

Product Datasheet

**Anti-Phospho-Insulin Receptor Beta (Tyr1185)/APC
[Out of stock] (orb1007360)**

Description

Anti-Phospho-Insulin Receptor Beta

Species/Host

Rabbit

Conjugation

APC

Tested

IF

Applications
Immunogen

KLH conjugated Synthesised phosphopeptide derived from human INSR around the phosphorylation site of Tyr1185 [DI(p-Y)ET]

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

68/152kDa

Uniprot ID
P06213, P15208, P15127
Dilution Range

IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

Expiration Date

12 months from date of receipt.