



Product Datasheet

PRKAB1 antibody (orb10077)

Description	PRKAB1 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human AMPK beta 1 (201-270/270 aa)
Target	AMPKB1 + PRKAB2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	30 kDa
Uniprot ID	Q9Y478
Expiration Date	12 months from date of receipt.



IHC-P of mouse colon carcinoma (AMPK bet...



Immunohistochemical analysis of formalin...