

Product Datasheet

SA1 protein stromal antigen 1 antibody (orb100818)

Description	SA1 protein stromal antigen 1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Nuclear protein stromal antigen 1 (901-1000/1258 aa)
Target	Nuclear protein stromal antigen 1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	138 kDa
Uniprot ID	Q8WVM7
Expiration Date	12 months from date of receipt.