


Product Datasheet

Human secretory IgA antibody (orb100834)

Description	Human secretory IgA antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Purified human sIgA
Target	human secretory IgA
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P01876
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Sample:
Lane 1:
U2os
(Human)
Cell
Lysate...