


Product Datasheet

beta Amyloid antibody (orb10087)


Description	beta Amyloid antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	FC
Immunogen	KLH conjugated synthetic peptide of human beta-Amyloid (1-42/42 aa)
Target	beta Amyloid 1-42
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	4.4 kDa
Uniprot ID	P05067
Dilution Range	Flow-Cyt=1µg/Test
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



Western blot of mouse tissue lysate usin...



IHC-P of rat hippocampus tissue (orb1008...