

Product Datasheet

LOXL2 antibody (orb100924)

Description

LOXL2 antibody

Species/Host

Rabbit

Reactivity

Bovine, Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human LOXL2. Please contact us for the exact immunogen sequence. The peptide is available as orb374805.

Target

LOXL2

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

87 kDa

Uniprot ID

[B5DF27](#), [Q9Y4K0](#), [P58022](#)

NCBI

[21](#), [002318](#), [002309](#), [11](#)

Entrez


[4017](#)

Dilution Range


WB:1:5000, IHC-P:1:500, IF/ICC:1:500

Expiration Date


12 months from date of receipt.



WB analysis of rat uterus tissue using L...



Western blot analysis of human breast ca...



Immunohistochemical analysis of formalin...