

Product Datasheet

PLB (phospho-Ser16) antibody (orb100953)

Description	PLB (phospho-Ser16) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Cardiac phospholamban around the phosphorylation site of Ser16 (RA(p-S)TI)
Target	phospho-PLB (Ser16)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	5.7 kDa
Uniprot ID	P26678
Expiration Date	12 months from date of receipt.