



Product Datasheet

Cytokeratin 12 antibody (orb100977)

Description	Cytokeratin 12 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CK12/Cytokeratin 12 (151-250/494 aa)
Target	CK12
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	54 kDa
Uniprot ID	Q99456
Expiration Date	12 months from date of receipt.



Western blot analysis
of rat liver lysat...



Immunohistochemical
analysis of formalin...