


Product Datasheet

PDGFC antibody (orb101058)

Description	PDGFC antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human PDGFC/VEGF E (201-300/345 aa)
Target	PDGFC
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	39 kDa
Uniprot ID	Q9NRA1
Dilution Range	ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse stomach l...