

Product Datasheet

Myosin (H chain 1) antibody (orb101070)

Description	Myosin (H chain 1) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human Myosin heavy chain 1 (901-1000/1939 aa)
Target	Fast Myosin Skeletal Heavy chain
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	223 kDa
Uniprot ID	P12882
Expiration Date	12 months from date of receipt.