


Product Datasheet

Desmoglein 3 antibody (orb101121)

Description	Desmoglein 3 antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Desmoglein 3 (50-75/999 aa)
Target	Desmoglein 3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	105 kDa
Uniprot ID	P32926
Expiration Date	12 months from date of receipt.



Paraffin-embedded rat tongue tissue fixe...