



Product Datasheet

EphA8 antibody (orb101151)

Description	EphA8 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human EphA8/Eph receptor A8 (451-550/1005 aa)
Target	EphA8
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	108 kDa
Uniprot ID	P29322
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat brain...



Immunohistochemical staining of rat brain...