

Product Datasheet

BMP11 antibody (orb101175)

Description

BMP11 antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human BMP11. Please contact us for the exact immunogen sequence. The peptide is available as orb393924.

Target

BMP11

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

45 kDa

Uniprot ID

[Q9Z217](#), [O95390](#), [Q9Z1W4](#)

NCBI

[005811](#), [005802](#)

Entrez



[10220](#)

Dilution Range


IHC-P: 1:100-600, WB: 1:200-2000, WB: 1:200-1000

Expiration Date


12 months from date of receipt.



IHC-P staining of rat colon tissue using...



Immunohistochemical staining of paraffin...



Western blot analysis of mouse colon (la...