


Product Datasheet

FASTKD5 antibody (orb101178)

Description	FASTKD5 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from mouse FASTKD5 (501-600/762 aa)
Target	FASTKD5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	84 kDa
Uniprot ID	Q7TMV3
Expiration Date	12 months from date of receipt.



Immunohistochemical
staining of human
pr...