



Product Datasheet

GPR55 antibody (orb101191)


Description	GPR55 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human GPR55. Please contact us for the exact immunogen sequence. The peptide is available as orb374801.
Target	GPR55
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	37 kDa
Uniprot ID	F1MAK4, Q14BV9, Q9Y2T6
NCBI	21, 005674, 31, 005683
Entrez	9290
Dilution Range	WB: 1:100, IHC-P: 1:100, IF/ICC: 1:100
Expiration Date	12 months from date of receipt.



IHC-P image of mouse brain tissue using ...



Western blot analysis of Rat ovary tissue...



WB analysis of Rat kidney tissue using G...