

Product Datasheet

**Anti-Phospho-Dab1 (Tyr198) /Gold [Out of stock]
(orb1012377)**

Description

Anti-Phospho-Dab1 (Tyr198) /Gold

Species/Host

Rabbit

Conjugation

Gold

Tested

GICA, IEM, IGS

Applications
Immunogen

KLH conjugated synthesised phosphopeptide derived from human Dab1 around the phosphorylation site of Tyr198 [V(p-Y)QY]

Preservatives

0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

0.4mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

65kDa

Uniprot ID
O75553, P97318, Q8CJH2
Dilution Range

IEM=1:20-200 IGS=1:20-200
GICA=1:20-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.

Expiration Date

12 months from date of receipt.