


Product Datasheet

Syndecan 4 antibody (orb101265)

Description	Syndecan 4 antibody
Species/Host	Rabbit
Reactivity	Bovine, Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthetic peptide derived from human Syndecan 4 (101-198/198 aa)
Target	Syndecan 4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	20 kDa
Uniprot ID	P31431
Expiration Date	12 months from date of receipt.



Western blot analysis of rat brain lysat...