



Product Datasheet

STMN3 antibody (orb101319)


Description	STMN3 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human STMN3 (77-130/180 aa)
Target	STMN3
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21 kDa
Uniprot ID	Q9NZ72
Expiration Date	12 months from date of receipt.



Western blot analysis
in rat brain tissu...



Immunohistochemical
analysis of paraffin...



Immunohistochemical
analysis of paraffin...