



Product Datasheet

PP2A beta antibody (orb101330)


Description	PP2A beta antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human PP2A beta/PPP2CB (21-120/309 aa)
Target	PP2A alpha + beta
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36 kDa
Uniprot ID	P62714
Dilution Range	WB=1:500-2000, IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P analysis of Human colon cancer tissue using PP2A beta antibody.



IHC-P analysis of rat brain using PP2A beta antibody.



Western blot: Brain (Mouse); Primary: Anti-PP2A beta antibody.