



Product Datasheet

EIF2C1 antibody (orb101399)

Description	EIF2C1 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human EIF2C1/Ago1 (431-530/984 aa)
Target	EIF2C
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	94 kDa
Uniprot ID	Q9UL18
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse brain...



Immunohistochemical staining of human re...