



Product Datasheet

CLIP1 antibody (orb101438)

Description	CLIP1 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human CLIP1/CLIP170 (451-550/1438 aa)
Target	CLIP170
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	158 kDa
Uniprot ID	P30622
Expiration Date	12 months from date of receipt.



Flow cytometric analysis of HeLa cells
U...



Immunohistochemical staining of mouse
CO...