



Product Datasheet

ICAM1 antibody (orb101441)


Description	ICAM1 antibody
Species/Host	Rabbit
Reactivity	Goat, Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ICAM1. Please contact us for the exact immunogen sequence. The peptide is available as orb394009.
Target	ICAM1
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	58 kDa
Uniprot ID	P05362 , Q00238
Entrez	3383
Dilution Range	IHC-P: 1:50-300, IF/ICC: 1:50-300, WB: 1:200-2000
Expiration Date	12 months from date of receipt.



IF analysis of rat lymph node tissue using ICAM1 antibody.



IHC-P image of rat lymph node tissue using ICAM1 antibody.



Western blot analysis of mouse spleen (lane 1) using ICAM1 antibody.