


Product Datasheet

ASPP2 antibody (orb10146)

Description	ASPP2 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human ASPP2 (1001-1128/1128 aa)
Target	ASPP2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	124 kDa
Uniprot ID	Q13625
Dilution Range	IHC-P=1:100-500, ICC=1:50
Expiration Date	12 months from date of receipt.



Western Blot of Rat liver tissue using A...