


Product Datasheet

Angiotensinogen antibody (orb10147)

Description	Angiotensinogen antibody
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human Angiotensinogen (51-110/477 aa)
Target	Angiotensinogen
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49 kDa
CAS Number	3
Uniprot ID	P01019
Expiration Date	12 months from date of receipt.



Paraffin-embedded rat liver tissue fixed...