

Product Datasheet

**Anti-Beta galactosidase/PE-Cy5 [Out of stock]
(orb1014849)**

Description

Anti-Beta galactosidase/PE-Cy5

Species/Host

Rabbit

Reactivity

Human, Mouse, Rat

Conjugation

PE/Cy5

Tested

ICC

Applications
Immunogen

KLH conjugated synthetic peptide derived from human Beta galactosidase (221-290aa)

Target

GLB1

Preservatives

0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.

Form/Appearance

Lyophilized or Liquid

Concentration

1mg/ml

Storage

Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

MW

71kDa

Uniprot ID
P16278
Expiration Date

12 months from date of receipt.