



Product Datasheet

p73 alpha antibody (orb101487)

Description	p73 alpha antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from mouse p73 alpha (301-350/636 aa)
Target	P73
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70 kDa
Uniprot ID	Q9JJP2
Dilution Range	WB=1:500-2000, IHC-P=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
human
esophagus
using ...



IHC-P
analysis of
rat skin using
p73 alp...