



Product Datasheet

ENaC Alpha antibody (orb101489)


Description	ENaC Alpha antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human ENaC Alpha. Please contact us for the exact immunogen sequence. The peptide is available as orb12930.
Target	ENaC Alpha
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	76 kDa
Uniprot ID	Q61180 , P37089
NCBI	001029
Dilution Range	WB: 1:100-800, IHC-P: 1:100-600, IF/ICC: 1:100-600
Expiration Date	12 months from date of receipt.



IHC-P
staining of
rat pancreas
tissue us...



IF image of
guinea pig
kidney tissue
usi...



IHC-P
staining of
rat lung
tissue using
...