


Product Datasheet

SLC26A4 antibody (orb101538)

Description	SLC26A4 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human Solute carrier family 26 member 4 (301-400/780 aa)
Target	SLC26A4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	93 kDa
Uniprot ID	O43511
Dilution Range	WB=1:500-2000, Flow-Cyt=1ug/test
Expiration Date	12 months from date of receipt.



IHC-P of
human
colorectal
cancer using
S...