

Product Datasheet

Anti-Phospho-NMDAR2B (Tyr1474)/RBITC [Out of stock] (orb1015804)

Description

Anti-Phospho-NMDAR2B

Species/Host	Rabbit
Conjugation	RBITC
Tested Applications	ICC, IF
Immunogen	KLH conjugated synthesised phosphopeptide derived from human NMDAR2B around the phosphorylation site of Tyr1474
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	1mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	163kDa
Uniprot ID	Q00960, Q01097, Q13224
Dilution Range	ICC=1:50-200 IF=1:50-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.