

Product Datasheet

AVPR2 antibody (orb10161)

Description	AVPR2 antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	IF, IHC-Fr, WB
Immunogen	KLH conjugated synthetic peptide derived from human AVPR2 (281-371/371 aa)
Target	AVPR2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40 kDa
Uniprot ID	P30518
Expiration Date	12 months from date of receipt.