

Product Datasheet

NIFK (phospho-Thr234) antibody (orb101615)

Description	NIFK (phospho-Thr234) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthesised phosphopeptide derived from human NIFK/MKI67IP around the phosphorylation site of Thr234 (GP(p-T)PV)
Target	phospho-NIFK (Thr234)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34 kDa
Uniprot ID	Q8TB66
Expiration Date	12 months from date of receipt.