



Product Datasheet

CREB1 antibody (orb10174)

Description	CREB1 antibody
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human CREB-1 (101-150/341 aa)
Target	CREB-1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	37 kDa
Uniprot ID	P16220
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



IHC-p of the mouse hippocampus shows CRE...



A single band of CREB-1 is detected in m...