

Product Datasheet

GATA-4 (phospho-Ser262) antibody (orb101765)

Description	GATA-4 (phospho-Ser262) antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human GATA4 around the phosphorylation site of Ser262 (RL(p-S)AS)
Target	phospho-GATA4 (Ser262)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49 kDa
Uniprot ID	P43694
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.