

Product Datasheet

Viral ODC1 protein (orb10178)

Description	Recombinant of viral ODC1 protein
Endotoxins	Not test.
Preservatives	If the delivery form is liquid, the default storage buffer is Tris/PBS-based buffer, 5%-50% glycerol. If the delivery form is lyophilized powder, the buffer before lyophilization is Tris/PBS-based buffer, 6% Trehalose, pH 8.0.
Form/Appearance	Liquid or Lyophilized powder
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Tag	N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged
Note	For research use only
Application notes	Tag info: N-terminal 10xHis-SUMO-tagged and C-terminal Myc-tagged Expression Region: 1-461aa Sequence Info: Full Length Glycerol content: 0.5
Protein Sequence	MNNFGNEEFDCFLDEGFTAKDILDQKINEVSSDDKDAFYVADLGDILKKHLR WLKALPRVTPFYAVKCNDSKAIVKTLAATGTGFDCASKTEIQLVQSLGVPPERIIY ANPCKQVSQIKYAANNGVQMMTFDSEVELMKVARAHPKAKLVLRATDDSKAV CRLSVKFGATLRTSRLLLERAKELNIDVVGVSFHVSGSGCTDPETFVQAISDARCV FDMGAEVGFSGMYLLDIGGGFPGSEDEVKLKFEETGVINPALDKYFPSDSGVRIAE PGRYYVASAFTLAVNIIAKKIVLKEQTGSDDDESESEQTFMYVNDGVYGSFNCIL YDHAHVKPLLQKRPKPDEKYYSSSIWGPTCDGLDRIVERCDLPEMHVGDWMLF ENMGAYTVAAASTFNGFQRPTIYYVMSGPAWQLMQQFQNPDPPEVEEQDAST LPVSCAWESGMKRHRACASASINV
Purity	Greater than 90% as determined by SDS-PAGE.
Source	E.coli
MW	71.1 kDa
Uniprot ID	P11926
Expiration Date	6 months from date of receipt.