

Product Datasheet

Beclin 1 antibody (orb10182)

Description

Beclin 1 antibody

Species/Host

Rabbit

Reactivity

Human, Rat

Conjugation

Unconjugated

Tested

ICC, IF, IHC-P, WB

Applications

Immunogen

KLH conjugated synthetic peptide derived from human Beclin 1. Please contact us for the exact immunogen sequence. The peptide is available as orb12882.

Target

Beclin 1

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

52 kDa

Uniprot ID

[Q14457](#), [Q91XJ1](#), [O88597](#)

NCBI

[003766](#), [31](#), [003757](#), [11](#)

Entrez



[8678](#)

Dilution Range


WB: 1:100-2000, IHC-P: 1:100-600, IF/ICC: 1:100-600

Expiration Date


12 months from date of receipt.



Immunohistochemical staining of rat spin...



Immunohistochemical staining of paraffin...



IF image of rat lymph node tissue using ...