


Product Datasheet

Beta casein antibody (orb10184)

Description	Beta casein antibody
Species/Host	Rabbit
Reactivity	Bovine, Goat, Human, Sheep
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from sheep Beta-casein (151-222/222 aa)
Target	Beta casein
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	23 kDa
Uniprot ID	P11839
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of sheep non-fat m...