


Product Datasheet

Casein antibody (orb10185)

Description	Casein antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Purified casein protein was used as immunogen to make this antibody
Target	Casein
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22 kDa
Dilution Range	WB: 1:100-1000 (based on 0.5mg/ml)
Expiration Date	12 months from date of receipt.



Western blot analysis of bovine milk (la...