



Product Datasheet

ZNF185 antibody (orb101869)

Description	ZNF185 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from human ZNF185 (551-650/689 aa)
Target	ZNF185
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	76(H)/39(M) kDa
Uniprot ID	O15231
Dilution Range	WB=1:500-2000, Flow-Cyt=3ug/Test
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Rat panc...



Western blot analysis of extracts from P...