



Product Datasheet

LCMT1 antibody (orb101882)


Description	LCMT1 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human LCMT1 (45-150/334 aa)
Target	LCMT1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	38 kDa
Uniprot ID	Q9UIC8
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse Pancreas ...



Western blot analysis of mouse Cerebellu...



Immunohistochemical staining of rat test...