

Product Datasheet

MCK10 (phospho-Tyr513) antibody (orb101885)

Description	MCK10 (phospho-Tyr513) antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human MCK10 around the phosphorylation site of Tyr513 (PA(p-Y)RL)
Target	phospho-MCK10 (Tyr513)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	99 kDa
Uniprot ID	Q08345
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.