

Product Datasheet

Anti-phospho-PTK9 (Tyr321)/Gold [Out of stock] (orb1018888)

Description	Rabbit polyclonal antibody to Anti-phospho-PTK9 (Tyr321)/Gold
Species/Host	Rabbit
Conjugation	Gold
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human TWF1 around the phosphorylation site of Tyr321
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	40kDa
Uniprot ID	Q56JV6 , Q12792 , Q91YR1 , Q5RJR2
Dilution Range	IEM=1:20-200 ICA=1:20-200 ChIP=1:20-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.