

Product Datasheet

Anti-phospho-MEK1 + MEK2 (Ser222 + Ser226)/Gold [Out of stock] (orb1018987)

Description	Rabbit polyclonal antibody to Anti-phospho-MEK1 + MEK2 (Ser222 + Ser226)/Gold
Species/Host	Rabbit
Conjugation	Gold
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human MEK1 around the phosphorylation site of Ser222
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43kDa
Uniprot ID	Q02750, P36507
Dilution Range	IEM=1:20-200 ICA=1:20-200 ChIP=1:20-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.