


Product Datasheet

BIM antibody (orb10190)

Description	BIM antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	KLH conjugated synthetic peptide derived from the middle of human Bim (51-150/198 aa)
Target	Bim
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	22 kDa
Uniprot ID	043521
Dilution Range	WB=1:500-2000, Flow-Cyt=1ug/test
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...