


Product Datasheet

BMPR1B antibody (orb101902)

Description	BMPR1B antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human BMPR1B (61-160/502 aa)
Target	BMPR1B
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	56 kDa
Uniprot ID	000238
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.



Western blot:U251(human)Cell;U87mg(human...)