

Product Datasheet

Anti-phospho-EGFR (Tyr1092)/Gold [Out of stock] (orb1019042)

Description	Anti-phospho-EGFR (Tyr1092)/Gold
Species/Host	Rabbit
Conjugation	Gold
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human EGFR around the phosphorylation site of Tyr1092
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	170kDa
Uniprot ID	P00533, Q01279
Expiration Date	12 months from date of receipt.