

Product Datasheet

Anti-phospho-Tau (Thr231)/Gold [Out of stock] (orb1019167)

Description	Anti-phospho-Tau (Thr231)/Gold
Species/Host	Rabbit
Conjugation	Gold
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Tau protein isoform 2 around the phosphorylation site of Thr231 [VR(p-T)PP]
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52/79kDa
Uniprot ID	P10637, P29172, P10636, P19332
Expiration Date	12 months from date of receipt.