

Product Datasheet

**Anti-Phospho-p70 S6 Kinase Beta (Thr228)/Gold
[Out of stock] (orb1019406)**

Description

Anti-Phospho-p70 S6 Kinase Beta

Species/Host	Rabbit
Conjugation	Gold
Tested Applications	GICA, IEM, IGS
Immunogen	KLH conjugated synthesised phosphopeptide derived from human p70 S6 Kinase Beta around the phosphorylation site of Thr228
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	59kDa
Uniprot ID	Q9Z1M4, Q9UBS0
Dilution Range	IEM=1:20-200 IGS=1:20-200 GICA=1:20-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.