

Product Datasheet

Anti-Phospho-p70 S6 Kinase Beta 2 (Tyr389)/Gold [Out of stock] (orb1019407)

Description	Rabbit polyclonal antibody to Anti-Phospho-p70 S6 Kinase Beta 2 (Tyr389)/Gold
Species/Host	Rabbit
Conjugation	Gold
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human p70 S6 Kinase Beta 2 around the phosphorylation site of Tyr389
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at 2-8°C for 3-6 months. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53kDa
Uniprot ID	Q9Z1M4, Q9UBS0
Dilution Range	IEM=1:20-200 ICA=1:20-200 ChIP=1:20-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.