

Product Datasheet

Anti-Phospho-Insulin Receptor Beta (Tyr1355)/Gold [Out of stock] (orb1019413)

Description

Anti-Phospho-Insulin Receptor Beta

Species/Host	Rabbit
Conjugation	Gold
Tested Applications	GICA, IEM, IGS
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Insulin Receptor alpha around the phosphorylation site of Tyr1355 [RS(p-Y)EE]5
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68/152kDa
Uniprot ID	P06213, P15208, P15127
Dilution Range	IEM=1:20-200 IGS=1:20-200 GICA=1:20-200 not yet tested in other applications. optimal dilutions/concentrations should be determined by the end user.
Expiration Date	12 months from date of receipt.