

Product Datasheet

Anti-phospho-FAK(Ser732)/Gold (orb1020034)

Custom Reagents.

Anti-phospho-FAK(Ser732)/Gold

Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Gold
Tested Applications	ICC
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human FAK around the phosphorylation site of Tyr732
Target	PTK2
Preservatives	0.02M TBS(pH8.2) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Lyophilized or Liquid
Concentration	0.4mg/ml
Storage	Store at -20 °C for one year. Avoid repeated freeze/thaw cycles. The lyophilized antibody is stable at room temperature for at least one month and for greater than a year when kept at -20°C. When reconstituted in sterile pH 7.4 0.01M PBS or diluent of antibody the antibody is stable for at least two weeks at 2-4 °C.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	116kDa
Uniprot ID	Q05397
Expiration Date	12 months from date of receipt.